Venus® Color

DE	Gebruiksaanweisung	SE	Bruksanvisning
GB	Instructions for use	DK	Brugervejledning
FR	Mode d'emploi	NO	Bruksanvisning
ES	Instrucciones de uso	FI	Käyttöohjeet
IT	Istruzioni per l'uso		
PT	Instruções de uso		
NL	Gebruiksaanwijzing		

60654712/00

Venus® Color

Lichtbehandelend composiet voor individuele karakterisering en tot afdekking van Verfärbungen der Zahnoberflächen bei direktem und indirektem Komposit-Restaurationen. Farben: WHITE, BLUE, CORN, AMBER, CHOCO.

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Venus® Color

Light-curing composite for customized characterisation and covering of discolourations on hard tooth tissue in direct and indirect composite restorations. Colors: WHITE, BLUE, CORN, AMBER, CHOCO.

Composition
Venus Color is based on multifunctional methacrylate monomers (UDMA, EBADMA) and contains approximately 65% by weight or 38% by volume inorganic fillers, such as Ba-Al-F-silicate glass, YbF3, SiO2. The filler particle size is between 0.02 µm and 5 µm.

Areas of application
Customised characterisation of direct and indirect composite restorations, veneers, acrylic teeth in partial and complete prostheses and temporary appliances: enamel cracks, chally spots, abrasion surfaces, fissures, discolouration of tooth neck, enamel bleaches, coloured intermediate layer for covering discolourations.

Contraindications
Do not use in case of a known allergy to any of the constituents of Venus Color. Direct application on dentine of vital teeth without dentine protection / dentine bonding agent. If dry dentin or the prescribed application technique is not possible.

Side effects
This product or one of its components may in particular cause hypersensitive reactions. If suspected, information on the ingredients can be obtained from the manufacturer.

- Avoid contact with skin. If product comes into contact with the skin, wash thoroughly with soap and water immediately. May cause sensitisation by skin contact. Avoid contact with eyes. Rinse immediately and thoroughly with plenty of water in case of contact with eyes and contact an ophthalmologist.
- Do not allow it to come into contact with preparations containing eugenol or phenolic disinfectants, as they may restrict the optimal polymerisation of the stains.
- Visually inspect for damage before first use. Damaged products must not be used.
- The cannulas supplied are dosing aids and not intended for direct application on the patient.

Kind of application
Venus Color is employed in combination with composite materials or acrylic teeth (e.g. Pal® prosthetic teeth from Heraeus Kulzer). The ready-to-use Venus Color paints are applied undiluted in a thin layer or diluted with a transparent, flowable composite (e.g. Venus® Diamond Flow CL) in different techniques. Depending on the application, they can be applied indirectly with a brush or other suitable instruments (e.g. probe). To reduce the color intensity of the Venus Color stains, they can be applied in a thin layer or mixed with a transparent, flowable composite (e.g. Venus® Diamond Flow CL).

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